

Supplemental Table 1. Cut Slide Age (antigen stability) for ALK

Supplemental Table 1. ALK (D5F3) Cut Slide Age Scores for samples stored at 30°C											
Lot	Slide ID	Fresh Cut		Month 3.5		Month 4		Month 5		Month 6	
		Stain	Bkgd	Stain	Bkgd	Stain	Bkgd	Stain	Bkgd	Stain	Bkgd
#1	Tissue 1	SP (3)	0	SP (2)	0	SP (2)	0	SP (3)	0	SP (3)	0
	Tissue 2	MP (2.5)	0	MP (2)	0	MP (2)	0	MP (2.5)	0	MP (2.5)	0
	Tissue 3	WP (1)	0	WP (1)	0	N (0.5)	0	WP (1)	0	N (<1)	0
	Tissue 4	N	0	N	0	N (0)	0	N	0	N	0
	Negative Rabbit Control	N (0)	0	N (0)	0	N (0)	0	N (0)	0	N (0)	0
#2	Tissue 1	SP (3)	0	SP (2)	0	SP (2)	0	SP (3)	0	SP (3)	0
	Tissue 2	MP (2.5)	0	MP (2.5)	0	MP (2)	0	MP (2.5)	0	MP (2.5)	0
	Tissue 3	WP (1)	0	WP (1)	0	N (0.5)	0	WP (1)	0	WP (1)	0
	Tissue 4	N	0	N	0	N (0)	0	N	0	N	0
	Negative Rabbit Control	N (0)	0	N (0)	0	N (0)	0	N (0)	0	N (0)	0
#3	Tissue 1	SP (3)	0	SP (2)	0	SP (2)	0	SP (3)	0	SP (3)	0
	Tissue 2	MP (2.5)	0	MP (2.5)	0	MP (2)	0	MP (2.5)	0	MP (2.5)	0
	Tissue 3	WP (1)	0	WP (1)	0	N (0.5)	0	WP (1)	0	N (<1)	0
	Tissue 4	N	0	N	0	N (0)	0	N	0	N	0
	Negative Rabbit Control	N (0)	0	N (0)	0	N (0)	0	N (0)	0	N (0)	0

SP= Strong , MP= Moderate, WP=Weak, N= negative, Bkgd= Background