


Standard titration


Slow titration

Fast titration


A. LTG monotherapy



B. With concomitant medicine inhibiting glucuronic acid conjugation. (ex. VPA)



C. With concomitant medicine promoting glucuronic acid conjugation. (ex. CBZ)



Supplemental Digital Content 1. Examples of the titration of LTG.

Figures illustrating dosage and increasing intervals of standard titration schedule (green), samples from various slow (blue) and fast (yellow) titration schedules. Examples of the titration schedules of LTG monotherapy is shown in (A). Since LTG inhibits (B) or promotes (C) glucuronic acid conjugation, the titration methods of LTG differs depending on the concomitant medicine. Abbreviations: VPA, valproate; CBZ, carbamazepine.